




Saturn 5TM

Laser Systems





Saturn 5™ Lasers

Fixed and Directional Laser Systems

*“The moveable
Saturn™ laser
means I can ablate
the zona exactly
where I want to,
without having to
move the embryo”*

Samantha Knight,
SPIRE London Fertility Clinic, UK

Precise

Sub-micron accuracy and unique computer controlled laser with guaranteed laser alignment

Easy to Use

Intuitive RI Viewer™ software with streamlined user interface. An optional programmable foot pedal controls software and laser functions

Safest Power

Lowest laser pulse times for minimal energy near critical cells. Exclusion Zone™ feature ensures cell safety

Curved Biopsy Mode

Biopsy Mode allows accurate laser drilling along a drawn straight or curved line

Rapid

Faster than ever directional laser increases functionality and decreases procedure times

Multi-Pulse Mode

Rapid fire options give more control over procedures



Saturn 5™ Lasers

Faster biopsies


The Saturn 5™ Biopsy Mode is helping you to significantly improve biopsying methods; potentially reducing procedure time, lowering the incidence of blastocyst collapse and eliminating the need to mechanically tear off cells¹.

Using Biopsy Mode, you can draw a straight or curved line along the sample, select the number and size of holes, then simply “fire”. The laser will ablate exactly along the chosen path, meaning you no longer need to move the holding pipette at all. It is that easy.


The Saturn 5™ Biopsy Mode is quickly proving to be an essential tool in the practice of biopsying across the world.



Biopsy Mode



Directional Laser Target




Blastocyst Biopsy

Totally committed to embryo safety

Multiple safety features reassure you that the Saturn™ laser is the safest laser on the market. To keep your embryos safe, features include the Exclusion Zone² which can be set according to distance or temperature.

By employing a higher laser power, the Saturn™ Laser applies less total energy to make a specified hole size, in comparison to lower power laser systems³.

Energy delivered vs. hole diameter³ using Saturn™ Laser



Biopsy Mode allows ablation along a drawn curved or straight line

Always spot-on

As the Saturn 5™ pilot laser travels down the same fibre optic path as the ablation laser, it guarantees consistent positioning. You can calibrate and verify hole size and firing position with sub-micron accuracy simply and quickly, leaving more time for your procedures.

As the pioneer of moveable lasers for ART, the Saturn 5 Active™ laser is spearheading a revolution in hatching and biopsy techniques in labs all over the world.

Intuitive software as standard

In addition to its intuitive use and safety assurances, the new Saturn 5™ features RI Viewer™ software. The software offers uncluttered full screen imaging from the microscope, digital magnification and a modern, clean user interface. RI Viewer™ also offers a recording function, line measurement tools which are visible on-screen and are printable, plus a built-in laser simulator for training and demonstration purposes. RI Viewer™ supports a range of languages.

Applications⁴

For **Blastocyst/Trophectoderm Biopsy**, the Saturn 5 Active™ is unparalleled in its ease of use. The directional laser allows the user to make multiple ablations across the trophectoderm cells without needing to move the blastocyst. This gives the user superb accuracy, safety and incredible speed.

¹Lloyd S, Doshi A, Harper J, Application note. A new method of biopsying TE cells using the latest Saturn 5 Active™ Laser System, offers several potential ways to improve your procedures - Available on request.

²Chatzimeletiou, K., Picton, H.M & Handyside, A.H., 2001. Use of non-contact, infrared laser for zona drilling of mouse embryos: assessment of immediate effects on blastomere viability. Reproductive biomedicine online, 2(3), p.178. Available at: <http://www.ncbi.nlm.nih.gov/pubmed/12537793>.

Faster Biopsies

Sub-micron Accuracy

Safest Lasers on the Market

Intuitive Software

For **Blastomere Biopsy, Polar Body Biopsy and Blastocyst Collapsing (for vitrification)**, the Saturn 5 Active™ directional laser allows the embryo to stay in the desired position and focus so that ablations can be made wherever required without additional manipulation. Our Biopsy Mode also allows safe multi-pulse drilling along a predetermined line.

For **Assisted Hatching**, the directional laser means that accurate ablations can be made without the need to hold the embryo, making it very quick and accurate, with no additional consumable costs.

Using Saturn 5 Active™, you will find these procedures are almost effortless and can be performed quickly and accurately. These procedures can require difficult embryo manipulation when performed using fixed lasers. Using a Saturn 5 Active™ has clear advantages to the welfare of the embryo.

Saturn 5™ Laser Systems are Class I laser products as defined by international laser safety standards. They are CE-marked and FDA cleared⁵.

³RI White Paper. A comparison of different power levels used by laser systems in the IVF laboratory – Available upon request.

⁴The applicability of procedures is dependent on the regulations of the country into which the device is sold.

⁵In the USA, FDA cleared for clinical use for Laser Assisted Hatching (LAH) only.



Saturn 5™ Lasers

Specifications

Microscope Compatibility

Nikon

TMD, Diaphot 200/300
TE200/300, TE2000, Ti

Leica


DMIRB, DMI3000B/4000B/6000B
DMIL

Olympus

IMT2, IX50/70,
IX51/71/81, IX53/73/83

Zeiss

Axiovert 40/100/200,
Axio Observer



<i>Pilot Laser</i>	630-650nm spot targeting solid state diode laser - red pilot beam guarantees the position of the invisible ablation laser
<i>Objective</i>	Custom designed objective for optimum laser transmission, crystal clear imaging and minimal laser pulse times. Tested and proven not to exhibit astigmatism with Saturn 5™ Laser Systems. 20x and 40x objectives are available
<i>Ablation Laser</i>	1480nm / 400mW solid state diode laser. Pulse length range 0.001-2.0ms / 1-2000 μ s Class 1 laser product (IEC 60825-1:2007) Tested and proven not to exhibit thermal lensing
<i>Laser Unit Dimensions (WxDxH)</i>	220mm x 180mm x 34mm
<i>CRi Oosight™ Compatibility</i>	Compatible with CRi Oosight™ and SpindleView™ systems
<i>Fluorescence Compatibility</i>	Saturn™ Laser Systems are compatible with epi-fluorescence on selected microscopes. Contact us for details
<i>Operation Software</i>	RI Viewer™ imaging software included – with digital laser targeting
<i>PC System Requirements</i>	Operating systems: Windows 8.1, Windows 8, Windows 7 (SP1), Windows Vista, Windows XP (SP3)
<i>Mains Input</i>	100-240VAC, 50-60Hz



Research Instruments Ltd, Bickland Industrial Park, Falmouth, Cornwall, TR11 4TA, UK
Tel: +44 (0) 1326 372 753 Fax: +44 (0) 1326 378 783 e-mail: sales@research-instruments.com visit: www.research-instruments.com